

Sliding aids



EasyGlide 5020



EasyGlide 5060

EasyGlide

For sitting transfers from wheelchairs to beds, toilets, car seats and many other seated positions. It is especially suitable for leveling out height differences or for bridging short distances. Available in two models and several sizes.

Art. No.	Product name	Description	Unit
5020	EasyGlide	60 x 33 cm/ 23.6" x 13", three-part	1 Pce
5050	EasyGlide	75 x 33 cm/ 29.5" x 13", three-part	1 Pce
5060	EasyGlide	60 x 25 cm/ 23.6" x 9.8"	1 Pce
5070	EasyGlide	45 x 20 cm/ 17.7" x 7.9"	1 Pce

Accessories

4014	Storage pocket	For 5020, 5050	1 Pce
------	----------------	----------------	-------



ReadyGlide

ReadyGlide is a sliding board that has been designed to simplify transfers in a seated position, for example, between bed and wheelchair or between wheelchair and toilet.

Art. No.	Product name	Description	Unit
15030	ReadyGlide	60 x 25 cm/ 23.6" x 9.8", polyethylene	1 Pce



EasyGlide 5000

EasyGlide 5010

EasyGlide

For supine transfers from one bed to another.

Art. No.	Product name	Description	Unit
5000	EasyGlide, extra long	120 x 60 cm/ 47.2" x 23.6"	1 Pce
5010	EasyGlide, extra long	170 x 59 cm/ 66.9" x 23.2"	1 Pce



EasyGlide 5042

EasyGlide

For supine transfers from one bed to another but also for positioning higher up in bed and for placing a lifting sling.

Art. No.	Product name	Description	Unit
5040	EasyGlide, oval	85 x 39 cm/ 33.5" x 15.4"	1 Pce
5042	EasyGlide, oval	85 x 39 cm/ 33.5" x 15.4"	2 Pcs

EasyGlide, oval is also available in a smaller model for positioning further back in chair, see page 9.

