



OneWaySlide 2075



OneWaySlide 2069

OneWaySlide

OneWaySlide has many functions. The main one improves sitting positions in chairs or wheelchairs. The sliding mat effectively retards forward movement yet slides backwards, making it easy for the caregiver to help the user move into a more comfortable position in the chair. Also available in a tube-shaped design without handles for the more active user.


Art. No.	Product name	Description	Unit
2074	OneWaySlide	110 x 30 cm / 43.3" x 11.8", 6 handles	1 Pce
2075	OneWaySlide	130 x 40 cm / 51.2" x 15.7", 6 handles	1 Pce
2076	OneWaySlide	130 x 45 cm / 51.2" x 17.7", 6 handles	1 Pce
2077	OneWaySlide	130 x 55 cm / 51.2" x 21.7", 6 handles	1 Pce
2069	OneWaySlide, tubular	50 x 40 cm / 19.7" x 15.7"	1 Pce
2070	OneWaySlide, tubular	50 x 45 cm / 19.7" x 17.7"	1 Pce



EasyGlide

EasyGlide are used in pairs to facilitate positioning further back in chair. The smaller size makes it easier to place the sliding boards under the user's thighs and seat.

Art. No.	Product name	Description	Unit
5045	EasyGlide, oval	45 x 20 cm / 17.7" x 7.9"	2 Pcs



AutoSlide

Facilitates getting in and out of cars, allowing its users to easier turn and move sideways. The seat cushions are firmly attached to the backs of car seats, so users cannot slide forward.

Art. No.	Product name	Description	Unit
8100	AutoSlide	46 x 46 cm / 18.1" x 18.1"	1 Pce

