

## Sit-to-stand aids



### QuickMove

QuickMove is a safe and easy-to-use aid for sit-to-stand and transfer of users with impaired balance and reduced lower limb strength. Based upon natural movement patterns, QuickMove supports the user when rising, using his/her own strength, and then transfers the user, either standing up or sitting down, from bed to chair or wheelchair or to and from the toilet. The user must have some strength in the arms and be able to weightbear to actively participate in the procedure.

BackSupport is an accessory which can be used as torso/side support for QuickMove.

BackSupport provides extra stability and security for users during transfers with QuickMove.

Another alternative is SupportStraps which, together with EasyBelt/FlexiBelt/FlexiMove, can be used as torso/side support.

Art. No.	Product name	Description	Unit
400801334	QuickMove	Max 170 kg/ 375 lbs	1 Pce

#### Accessories

70200037	BackSupport	Torso/side support for QuickMove	1 Pce
6012 -- 6015	EasyBelt	Extra support for QuickMove	1 Pce
6002 -- 6006	EasyBelt, hug	Extra support for QuickMove	1 Pce
6113 -- 6116	EasyBelt, wipeable	Extra support for QuickMove	1 Pce
6023 -- 6026	FlexiBelt	Extra support for QuickMove	1 Pce
6072 -- 6076	FlexiBelt, hug	Extra support for QuickMove	1 Pce
3012, 3020, 3030	FlexiMove	Extra support for QuickMove	1 Pce
6061	SupportStraps	Accessory for belts	1 Pair



### FlexiMove

FlexiMove has multiple areas of application. For example as support when helping someone get up from a bed or chair, when moving higher up in bed or when transferring to and from a wheelchair – in that case, often on a sliding board. FlexiMove is available in two variants with different inside material.

Art. No.	Product name	Description	Unit
3012	FlexiMove	20 x 60 cm/ 7.9" x 23.6", nylon	1 Pce
3020	FlexiMove	20 x 50 cm/ 7.9" x 19.7", polyester velvet	1 Pce
3030	FlexiMove	20 x 60 cm/ 7.9" x 23.6", polyester velvet	1 Pce

#### Accessories

6061	SupportStraps	Accessory for belts and FlexiMove together with QuickMove	1 Pair
------	---------------	---	--------

